



Product Brochure





360° Protection

As a healthcare professional, you know how important it is to keep yourself and your patients safe from pathogens. That's why Hygiene360's product portfolio consists of a comprehensive range of fast-acting disinfectants for dental clinics, hospitals, polyclinics, emergency and healthcare services, as well as medical centres.




Infection Control in Healthcare Facilities: An On-going Battle!

Nosocomial infections, also known as HAIs (Hospital-Acquired Infections or Healthcare-Associated Infections) are a major concern for the healthcare industry. The European Centre for Disease Prevention and Control (ECDC) estimated that one in 18 patients in the European hospitals developed healthcare-associated infection each day. Furthermore, based on the European Centre for Disease Prevention and Control report, there was an average prevalence of 7.1% in the European countries where 4,131,000 patients were affected by approximately 4,544,100 episodes of healthcare-associated infections every year.

Healthcare facilities have a higher risk of transmitting microorganisms, and, therefore, as a provider of highly effective infection control solutions, Hygiene360 offers a complete range of disinfectants and cleaning products for the said healthcare facilities.

Colour-Coded for Specific Applications Areas

When it comes to infection control and prevention, nothing must be left to chance. This is why Hygiene360's comprehensive, quick, and safe products has been colour-coded according to their application.



Product Overview

	Ready-to-use	Concentrate	200 ml	1 litre	2 litres	5 litres	14.5 x 20 cm (classic scent)	14.5 x 20 cm (lemon scent)	14.5 x 40 cm (classic scent)	19.3 x 21 cm (lemon scent)
Surfaces										
PROSEPT® Spray	X			X		X				
PROSEPT® Ultima	X			X	X	X				
PROSEPT® Sensi	X		X	X		X				
PROSEPT® Wipes	X						120 wipes	120 wipes		
PROSEPT® Wipes Jumbo	X								100 wipes	
PROSEPT® Wipes Value	X							250 wipes		
PROSEPT® Wipes Premium	X									120 wipes
PROSEPT® Floor		X			X	X				
Instruments										
PROSEPT® Fortis		X		X	X	X				
PROSEPT® Instru		X		X	X	X				
PROSEPT® Burs	X			X	X	X				
PROSEPT® Sonic		X		X						
Dental Impressions and Suction Units										
PROSEPT® Jet Forte		X			X	X				
PROSEPT® Jet D		X			X	X				
PROSEPT® Impression S	X			X						

Product Overview

	Ready-to-use	100 ml	500 ml bottle with pump	500 ml (EURO)	1 litre (EURO)	5 litres
Hands						
PROSEPT® Gel	X	X	X	X	X	X

	1.3 ml output	150 ml	1 litre	3 litres	5 litres	10 litres	14.5 x 20 cm (120 pcs)	14.5 x 20 cm (250 pcs)	14.5 x 40 cm (100 pcs)	19.3 x 21 cm (120 pcs)
Accessories and Application Aids										
PROSEPT® Burs Container		X								
PROSEPT® Soaking Tray			X	X						
Trigger Sprayer I1 (1.3 ml)	X		X							
Canister Tap					X	X				
PROSEPT® Wipes T1 (empty tub)							X			
PROSEPT® Wipes T2 (empty tub)								X	X	
PROSEPT® Wipes T3 (empty tub)										X



**Ready-to-Use Disinfectants and
Cleaners for Small Surfaces**

PROSEPT® Spray

Ready-to-use solution for the quick and comprehensive disinfection of non-invasive medical devices.

- **Alcohol content of less than 50%**

The disinfectant is particularly suitable for use on sensitive surfaces such as artificial leather covers.

- **Ergonomically designed bottle**

The ergonomically designed bottle ensures a safe and easy application of the disinfectant.

Instructions for use: For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds, auxiliaries and water.

Delivery forms: 1 litre bottle, 5 litre canister.

Complies with EN 13624, EN 13727, EN 14348, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	30 seconds
Poliovirus	60 seconds



PROSEPT® Ultima

Ready-to-use solution for the quick disinfection of non-invasive medical devices.

- **Aldehyde-, chlorine- and phenol-free**

Comprehensive disinfection with no health risks.

- **Alcohol content of less than 23%**

The disinfectant is particularly suitable for use on sensitive surfaces such as artificial leather covers or acrylic glass.

Instructions for use: For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 7.5 g ethanol, 15 g isopropanol, 0.76 g quaternary ammonium compounds, auxiliaries and water.

Delivery forms: 1 litre bottle, 2 litre bottle, 5 litre canister.

Complies with EN 13624, EN 13727, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Limited spectrum virucidal	5 minutes
Adenovirus	5 minutes
Norovirus	5 minutes



PROSEPT® Sensi

Alcohol-free, ready-to-use solution for the quick and comprehensive disinfection of non-invasive medical devices.

- **Alcohol-free**

The alcohol-free disinfectant is suitable for the disinfection of sensitive surfaces and objects such as acrylic glass.

- **Long-lasting antimicrobial action**

The antimicrobial action remains after the product has dried, thus, effectively reducing the risk of recontamination.

Instructions for use: For your protection, wear disposable gloves and safety goggles. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. Wipe medical devices, which come into contact with skin, once a day with a water-soaked single-use cloth to prevent the build-up of biocides. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 0.4 g benzalkonium chloride, 0.4 g N-alkyl aminopropyl glycine, auxiliaries and water.

Delivery forms: 200 ml foamer bottle, 1 litre bottle, 5 litre canister.

Complies with EN 13624, EN 13727, EN 14348, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 minutes
Mycobactericidal	60 minutes
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 minutes
Adenovirus	30 minutes
Norovirus	60 minutes
Poliovirus	60 minutes
Rotavirus	30 seconds



A teal floor buffer with a white microfiber pad is shown cleaning a wooden floor. The buffer is positioned at the top of the frame, and the pad is in contact with the floor. The floor has a natural wood grain. The image is framed by a white curved border on the left and bottom. The background is a solid green color.

**Disinfectant and Cleaning
Concentrates for Large
Surfaces**

PROSEPT® Floor

Concentrated liquid for the disinfection and cleaning of non-invasive medical devices.

- **Extremely economical concentrate**

The extremely economical concentrate significantly reduces your disinfection expenses.

- **Good material compatibility**

The product is suitable for the disinfection of a wide variety of water-washable surfaces.

- **Simple application**

Volumes are indicated on the lid for fast and easy dosing of disinfectant concentrates.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Completely wipe or immerse the surface of the medical device to be disinfected and leave to act. Adhere to contact time. If necessary, dry the surface with a clean single-use cloth. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 10 g quaternary ammonium compounds, auxiliaries and water.

Delivery forms: 2 litre bottle, 5 litre canister.

Complies with EN 13624, EN 13727, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	2% in 1 min
Antibiotic-resistant bacteria	2% in 1 min
Yeasticidal	2% in 1 min
Antibiotic-resistant yeasts	2% in 1 min
Virucidal against enveloped viruses	2% in 5 min
Bovine viral diarrhea virus (BVDV)	2% in 5 min
Coronavirus	2% in 5 min
Hepatitis B virus	2% in 5 min
Hepatitis C virus	2% in 5 min
Herpes simplex virus	2% in 5 min
Human immunodeficiency virus (HIV)	2% in 5 min
Influenza A virus	2% in 5 min
Vaccinia virus	2% in 5 min
Limited spectrum virucidal	2% in 5 min
Adenovirus	2% in 1 min
Norovirus	2% in 5 min



Dosing table

Cap (1 Cap = 20 ml)	1x ☐ (= 20 ml) 2%
Water	980 ml
Time	5 min





Ready-to-Use
Disinfectant Wipes

PROSEPT® Wipes

Ready-to-use wipes for the quick, comprehensive, and residue-free disinfection of non-invasive medical devices.

- **Ready-to-use disinfectant wipes**

The use of wipes eliminates the need for spray disinfection, which can cause respiratory tract, eye, and nose irritation.

- **Residue-free drying**

The disinfectant wipes leave no streaks and, therefore, make additional wiping redundant.

- **Dermatologically tested**

Instructions for use: For your protection, wear disposable gloves and safety goggles. Before use, cut open the refill pouch. Pull the first wipe at the heart of the roll through the star-shaped opening in the lid. Close the wipes box firmly after each use. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.


Composition in 100 g moisturizing solution: 31.3 g ethanol, 16.7 g isopropanol, auxiliaries and water.

Delivery forms: Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), classic scent, Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), lemon scent, PROSEPT® Wipes Tub T1 (empty tub).

Complies with EN 13624, EN 13727, EN 14348, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	60 seconds
Poliovirus	60 seconds
Rotavirus	30 seconds



PROSEPT® Wipes Jumbo

Ready-to-use wipes for the quick, comprehensive, and residue-free disinfection of non-invasive medical devices.

- **Ready-to-use wet wipes**

The use of wipes eliminates the need for spray disinfection, which can cause respiratory tract, eye, and nose irritation.

- **Residue-free drying**

The disinfectant wipes leave no streaks and, therefore, make additional wiping redundant.

- **Dermatologically tested**

Instructions for use: For your protection, wear disposable gloves and safety goggles. Before use, cut open the refill pouch. Pull the first wipe at the heart of the roll through the star-shaped opening in the lid. Close the wipes box firmly after each use. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g moisturizing solution: 31.3 g ethanol, 16.7 g isopropanol, auxiliaries and water.

Delivery forms: Refill pouch with 100 disinfectant wipes (14.5 x 40 cm), classic scent, PROSEPT® Wipes Tub T2 (empty tub).

Complies with EN 13624, EN 13727, EN 14348, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	60 seconds
Poliovirus	60 seconds
Rotavirus	30 seconds



PROSEPT® Wipes Value

Ready-to-use wipes for the quick, comprehensive, and residue-free disinfection of non-invasive medical devices.

- **Ready-to-use wet wipes**

The use of wipes eliminates the need for spray disinfection, which can cause respiratory tract, eye, and nose irritation.

- **Residue-free drying**

The disinfectant wipes leave no streaks and, therefore, make additional wiping redundant.

- **Dermatologically tested**

Instructions for use: For your protection, wear disposable gloves and safety goggles. Before use, cut open the refill pouch. Pull the first wipe at the heart of the roll through the star-shaped opening in the lid. Close the wipes box firmly after each use. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g moisturizing solution: 31.3 g ethanol, 16.7 g isopropanol, auxiliaries and water.

Delivery forms: Refill pouch with 250 disinfectant wipes (14.5 x 20 cm), lemon scent, PROSEPT® Wipes Tub T2 (empty tub).

Complies with EN 13624, EN 13727, EN 14348, EN 14476, and EN 16615.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	60 seconds
Poliovirus	60 seconds
Rotavirus	30 seconds



PROSEPT® Wipes Premium

Ready-to-use wipes for the quick, comprehensive, and residue-free disinfection of non-invasive medical devices.

- **Ready-to-use wet wipes**

The use of wipes eliminates the need for spray disinfection, which can cause respiratory tract, eye, and nose irritation.

- **Residue-free drying**

The disinfectant wipes leave no streaks and, therefore, make additional wiping redundant.

- **Dermatologically tested**

Instructions for use: For your protection, wear disposable gloves and safety goggles. Before use, cut open the refill pouch. Pull the first wipe at the heart of the roll through the star-shaped opening in the lid. Close the wipes box firmly after each use. Completely moisten and wipe the entire surface of the medical device. Leave disinfectant to act and adhere to contact time. Use on alcohol compatible materials only. Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g moisturizing solution: 31.3 g ethanol, 16.7 g isopropanol, auxiliaries and water.

Delivery forms: Refill pouch with 120 disinfectant wipes (19.3 x 21 cm), lemon scent PROSEPT® Wipes Tub T3 (empty tub)

Complies with EN 13624, EN 13727, EN 14348, EN 14476, and EN 16615.



Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	60 seconds
Poliovirus	60 seconds
Rotavirus	30 seconds



PROSEPT® Wipes

Effective Against 101 Pathogens

Effective in
60
Seconds



Compatible with 275 materials including artificial leather



Silicon star-shaped opening prevents wipes from drying out.



Durable lid



PROSEPT® Wipes
14.5 cm x 20 cm, 120 pcs



PROSEPT® Wipes Jumbo
14.5 cm x 40 cm, 100 pcs



















PROSEPT® Wipes Value
14.5 cm x 20 cm, 250 pcs



PROSEPT® Wipes Premium
19.3 cm x 21 cm, 120 pcs

PROSEPT® Wipes are compatible with more than 275 materials!



Original	After 2 weeks (Water)	After 2 weeks (PROSEPT® Wipes)	Description
			Material: Synthetic leather Tradename: KAVO® Time: 2 weeks reprocessing cycles Result: No change. Comparable to water.
			Material: Synthetic leather Tradename: Planmeca® Time: 2 weeks reprocessing cycles Result: No change. Comparable to water.
			Material: Synthetic leather Tradename: Castellini® Time: 2 weeks reprocessing cycles Result: No change. Comparable to water.
			Material: Synthetic leather Tradename: Anthos® Time: 2 weeks reprocessing cycles Result: No change. Comparable to water.
			Material: Synthetic leather Tradename: Joerg & Sohn Time: 2 weeks reprocessing cycles Result: No change. Comparable to water.



Want to find out more about our tests?
Call us at +41 (0)55 511 03 33 or write us an email at info@udm.swiss





Disinfection and Cleaning of Dental and Medical Instruments

PROSEPT® Fortis

Powerful concentrate for the disinfection and cleaning of medical instruments through manual reprocessing.

- **Contains effective corrosion inhibitors**

The disinfectant contains highly effective corrosion inhibitors to guard against rusting and it is less corrosive than water on some metals.

- **No discoloration or fixation of blood**

The solution does not cause discoloration or protein fixation on instruments.

- **Suitable for use in ultrasonic cleaning devices**

The cleansing power of the product can be further increased using an ultrasonic cleaning device.

- **Fully virucidal**

PROSEPT® Fortis is active against bacteria, mycobacteria, yeasts as well as enveloped and non-enveloped viruses.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Immerse instruments completely and adhere to contact time. After disinfection, rinse instruments thoroughly under running water of at least drinking quality for 10 seconds. Visually check instruments after the disinfection for contaminants. If necessary, remove contaminants with a brush and repeat the disinfection process. Replace working solution daily or when visibly contaminated. The disinfectant solution is also suitable for use in an ultrasonic cleaning device. Adhere to reprocessing recommendations of the instrument manufacturer. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 9.1 g alkylamine, 3.6 g didecyltrimethylammonium chloride, corrosion inhibitors, auxiliaries and water.

Delivery forms: 1 litre bottle, 2 litre bottle, 5 litre canister.


Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, and EN 17111.



Effective Spectrum

Bactericidal	3% in 5 min
Antibiotic-resistant bacteria	3% in 5 min
Tuberculocidal	3% in 5 min
Mycobactericidal	3% in 5 min
Yeasticidal	3% in 5 min
Antibiotic-resistant yeasts	3% in 5 min
Virucidal against enveloped viruses	3% in 5 min
Bovine viral diarrhea virus (BVDV)	3% in 5 min
Coronavirus	3% in 5 min
Hepatitis B virus	3% in 5 min
Hepatitis C virus	3% in 5 min
Herpes simplex virus	3% in 5 min
Human immunodeficiency virus (HIV)	3% in 5 min
Influenza A virus	3% in 5 min
Vaccinia virus	3% in 5 min
Virucidal	3% in 15 min
Adenovirus	3% in 10 min
Norovirus	3% in 10 min
Poliovirus	3% in 15 min

Dosing table

Cap (1 Cap = 20 ml)	1½x  (= 30 ml) 3%
Water	970 ml
Time	15 min



PROSEPT® Instru

Economical concentrate for the disinfection and cleaning of medical instruments through manual reprocessing.

- **Excellent material compatibility**

The disinfectant is compatible with a wide range of materials including aluminium.

- **Maximum cost-effectiveness**

The disinfectant is available in 5-litre canisters to reduce waste and offer better value for money.

- **Broad effective spectrum**

Active against all bacteria including those causing tuberculosis, yeasts, and viruses.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Immerse instruments completely and adhere to contact time. After disinfection, rinse instruments thoroughly under running water of at least drinking quality for 10 seconds. Visually check instruments after the disinfection for contaminants. If necessary, remove contaminants with a brush and repeat the disinfection process. Replace working solution daily or when visibly contaminated. The disinfectant solution is also suitable for use in an ultrasonic cleaning device. Adhere to reprocessing recommendations of the instrument manufacturer. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 4.55 g alkylamine, 1.8 g didecyltrimethylammonium chloride, corrosion inhibitors, auxiliaries and water.

Delivery forms: 1 litre bottle, 2 litre bottle, 5 litre canister.


Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, and EN 17111.



Effective Spectrum

Bactericidal	3% in 15 min
Antibiotic-resistant bacteria	3% in 15 min
Tuberculocidal	3% in 15 min
Mycobactericidal	3% in 15 min
Yeasticidal	3% in 15 min
Antibiotic-resistant yeasts	3% in 15 min
Virucidal against enveloped viruses	3% in 15 min
Bovine viral diarrhea virus (BVDV)	3% in 15 min
Coronavirus	3% in 15 min
Hepatitis B virus	3% in 15 min
Hepatitis C virus	3% in 15 min
Herpes simplex virus	3% in 15 min
Human immunodeficiency virus (HIV)	3% in 15 min
Influenza A virus	3% in 15 min
Vaccinia virus	3% in 15 min
Limited spectrum virucidal	3% in 60 min
Adenovirus	3% in 60 min
Norovirus	3% in 60 min

Dosing table

Cap (1 Cap = 20 ml)	1½x  (= 30 ml) 3%
Water	970 ml
Time	60 min



PROSEPT® Burs

Anti-corrosive, ready-to-use solution for the disinfection and cleaning of rotary dental and surgical instruments.

- **Fully virucidal in 30 seconds**

Active against even the most resistant viruses such as poliovirus and adenovirus.

- **Suitable for use in ultrasonic cleaning devices**

The cleansing power of the product can be further increased using an ultrasonic cleaning device.

- **Contains effective corrosion inhibitors**

The disinfectant contains highly effective corrosion inhibitors to guard against rusting and it is less corrosive than water on some metals.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Immerse the rotary instruments immediately after use in the disinfectant solution. Adhere to but do not exceed indicated contact time. After disinfection, leave the rotary instruments to dry. Do not rinse with water. Replace disinfectant solution daily or when visibly contaminated. Visually check instruments after the disinfection for contaminants. If necessary, remove contaminants with a brush and repeat the disinfection process. Replace working solution daily or when visibly contaminated. The disinfectant solution is also suitable for use in an ultrasonic cleaning device. Adhere to reprocessing recommendations of the instrument manufacturer. Use only with alkaline and alcohol compatible instruments. Not compatible with aluminium. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 32 g isopropanol, 0.7 g 2-phenoxyethanol, 0.2 g quaternary ammonium compounds, corrosion inhibitors, auxiliaries and water.

Delivery forms: 1 litre bottle, 2 litre bottle, 5 litre canister.

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, EN 14563, and EN 17111.



Effective Spectrum

Bactericidal	30 seconds
Antibiotic-resistant bacteria	30 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	30 seconds
Antibiotic-resistant yeasts	30 seconds
Fungicidal	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhoea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	30 seconds
Adenovirus	30 seconds
Norovirus	30 seconds
Poliovirus	30 seconds
Polyomavirus	30 seconds



PROSEPT[®] Sonic

Liquid concentrate for the thorough cleaning of contaminated dental and surgical instruments in the ultrasonic cleaning device.

- **Powerful cleanser**

The product's powerful cleansing properties make pre- and post-cleaning unnecessary.

- **Excellent cleansing properties**

Removes even stubborn residues such as blood, saliva, or plaster.

- **Simple application**

The bottle cap contains volume markings for fast and easy dosing.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Prepare the working solution according to the dosing table. Use an ultrasonic cleaning device for 30 minutes. After cleaning, rinse instruments thoroughly under running water of at least drinking quality for 10 seconds. Replace working solution daily or when visibly contaminated. Adhere to reprocessing recommendations of the instrument manufacturer. Not suitable for the cleaning of color-coded instruments and instruments that are fully or partly made of natural bristles, nylon, or aluminium. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 16 g potassium hydroxide, corrosion inhibitors, auxiliaries and water.

Delivery forms: 1 litre bottle.



Dosing table

Cap (1 Cap = 20 ml)	7½x ☐ (= 150 ml) 5%
Water	2850 ml
Time	30 min



Important Advice

- Do not use on colour-coded instruments and instruments made of/with natural bristles, nylon, or aluminium.
- Do not use metal brushes as they can damage the instrument alloy.
- The maximum standing time is 7 days. Depending on the degree of contamination, a daily change of the solution is recommended.





Disinfection of Dental Impressions

PROSEPT® Impression S

Ready-to-use spray disinfectant for dental impressions and impression trays.

- **Broad effective spectrum**

Active against all bacteria including those causing tuberculosis and viruses.

- **Short exposure time**

Active in 1 minute.

- **Will not alter the form or dimensional stability of dental impressions**

The disinfectant has no negative effect on the form or stability of dental impressions.

Instructions for use: For your protection, wear disposable gloves and safety goggles. Completely moisten the dental impression and impression tray and leave disinfectant to act. Adhere to but do not exceed indicated contact time. After disinfection, rinse impressions thoroughly under running water of at least drinking quality for 10 seconds. Use on alcohol compatible materials only. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds, auxiliaries and water.

Delivery forms: 1 litre bottle.

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 14561, EN 14562, and EN 14563.

Effective Spectrum

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	30 seconds
Poliovirus	60 seconds





Disinfection and Cleaning of Suction Units

PROSEPT® Jet D

Powerful liquid concentrate for the disinfection and cleaning of dental suction units.

- **Proven formula**

Dentists around the world are using this proven formula on a daily basis for the disinfection and cleansing of their suction units with complete satisfaction.

- **High-yielding concentrate**

2 litres of disinfectant concentrate yields 100 litres of working solution.

- **Simple application**

Volumes are indicated on the lid for fast and easy dosing of disinfectant concentrates.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Disinfect the suction unit daily or as per local regulations. Rinse the suction unit thoroughly with water after surgical procedures and prior to disinfection. Prepare 2 litres of working solution according to the dosing table. Aspirate 1.5 litres of the working solution and pour the remaining 500 ml down the spittoon bowl. Leave to act. Rinse thoroughly with water after disinfection. Not compatible with aluminium, copper, galvanized steel and polyurethane. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 2.475 g N-(3-aminopropyl)-N-dodecylpropane-1,3- diamine, 7.5 g quaternary ammonium compounds, anti-foaming agents, auxiliaries and water.

Delivery forms: 2 litre bottle, 5 litre canister.

Complies with EN 13624, EN 13727, EN 14476, EN 14561, and EN 14562.

Effective Spectrum

Bactericidal	2% in 15 min
Antibiotic-resistant bacteria	2% in 15 min
Yeasticidal	2% in 15 min
Antibiotic-resistant yeasts	2% in 15 min
Virucidal against enveloped viruses	2% in 15 min
Bovine viral diarrhoea virus (BVDV)	2% in 15 min
Coronavirus	2% in 15 min
Hepatitis B virus	2% in 15 min
Hepatitis C virus	2% in 15 min
Herpes simplex virus	2% in 15 min
Human immunodeficiency virus (HIV)	2% in 15 min
Influenza A virus	2% in 15 min
Vaccinia virus	2% in 15 min
Virucidal	2% in 60 min
Adenovirus	2% in 30 min
Norovirus	2% in 60 min
Poliovirus	2% in 30 min



Dosing table

Cap (1 Cap = 20 ml)	2x ☐ (= 40 ml) 2%
Water	1960 ml
Time	60 min



PROSEPT® Jet Forte

Liquid concentrate for the thorough cleaning of dental suction units.

- **High-yielding concentrate**
2 litres of concentrate yields 100 litres of working solution.
- **Excellent cleaning power**
Counteracts biofilm formation and bad odour, thereby, ensuring a reliable operation of the suction unit.
- **Non-foaming**
Gentle on dental suction units and hoses.
- **Simple application**
The bottle cap contains volume markings for fast and easy dosing.

Instructions for use: For your protection, wear a laboratory gown, household gloves and safety goggles. Clean the suction unit weekly or as per local regulations. Rinse the suction unit thoroughly with water after surgical procedures and prior to cleaning. Prepare 2 litres of working solution according to the dosing table. Aspirate 1.5 litres of the working solution and pour the remaining 500 ml down the spittoon bowl. Leave to act. Rinse thoroughly with water after cleaning. Not compatible with aluminium, copper, galvanized steel and polyurethane. A list of compatible and incompatible materials is available on our website. If you are uncertain as to the material composition of the medical device, contact the manufacturer before use. All serious accidents occurring in relation to the Medical Device should be reported to the manufacturer or the competent authority.

Composition in 100 g: 9.75 g alkaline cleaning compounds, anti-foaming agents, auxiliaries and water.

Delivery forms: 2 litre bottle, 5 litre canister.



Dosing table

Cap (1 Cap = 20 ml)	2x ☐ (= 40 ml) 2%
Water	1960 ml
Time	30 min





Hand Disinfectants

PROSEPT® Gel

Hydroalcoholic gel for the hygienic (30 seconds) and surgical (90 seconds) hand disinfection.

- **Broad effective spectrum**

Active against bacteria, mycobacteria, fungi, and viruses (including norovirus) in 60 seconds.

- **Non-sticky formulation**

The non-sticky formulation leaves no residue on your hands or gloves.

- **Contains anti-aging ingredients**

The hand disinfectant contains anti-aging ingredients, moisturizers, emollients, skin protection agents, vitamins, and soothing agents for maximum skin protection.

Instructions for use: Remove watches, rings, and other jewellery before applying the disinfectant. Apply required amount to dry and clean hands. For the hygienic hand disinfection: rub 3 ml for 30 seconds on hands and nails. For the surgical hand disinfection: rub twice 5 ml for 90 seconds on hands, nails, and forearms. Rub hands together in accordance with the 6 step technique for hand disinfection (EN 1500) and follow contact times.

Composition in 100 g: 52.8 g ethanol, 19.0 g 1-propanol, 2.8 g isopropanol, moisturizers, emollients, vitamin E acetate, panthenol, bisabolol, gelling agents.

Delivery forms: 100 ml bottle, 500 ml bottle, 500 ml bottle (EURO), 1 litre bottle (EURO), 5 litre canister.

Effective Spectrum

Bactericidal	15 seconds
Antibiotic-resistant bacteria	15 seconds
Tuberculocidal	15 seconds
Mycobactericidal	30 seconds
Yeasticidal	15 seconds
Antibiotic-resistant yeasts	15 seconds
Fungicidal	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	5 minutes
Adenovirus	60 seconds
Norovirus	15 seconds
Poliovirus	5 minutes
Rotavirus	15 seconds

Application area

Hygienic hand disinfection EN 1500	30 sec
Surgical hand disinfection EN 12791	90 sec



Easy Steps of Hands Disinfection



1 Apply PROSEPT® Gel to cupped hands.



2 Use the 6 step technique of hand disinfection pursuant to EN 1500. Palm to palm.



3 Right palm over the back of left hand (and vice versa).



4 Palm to palm with clasped fingers.



5 Top side of fingers against the palm of the opposite hand (both hands simultaneously).



6 Circular rubbing of the right thumb in the closed left fist (and vice versa).



7 Circular rubbing of the closed fingertips of the right hand in the left palm (and vice versa).



8 Repeat steps 2-7 for 30 seconds. The skin must remain moist during the entire exposure time. If necessary, coat the hands with PROSEPT® Gel again.



9 Observe the exposure time. 30 seconds for the hygienic hand disinfection and 90 seconds for surgical hand disinfection.



PROSEPT®

CONTENTS

- 1. Instrument
- 2. Instrument
- 3. Instrument
- 4. Instrument
- 5. Instrument
- 6. Instrument
- 7. Instrument
- 8. Instrument
- 9. Instrument
- 10. Instrument
- 11. Instrument
- 12. Instrument
- 13. Instrument
- 14. Instrument
- 15. Instrument
- 16. Instrument
- 17. Instrument
- 18. Instrument
- 19. Instrument
- 20. Instrument
- 21. Instrument
- 22. Instrument
- 23. Instrument
- 24. Instrument
- 25. Instrument
- 26. Instrument
- 27. Instrument
- 28. Instrument
- 29. Instrument
- 30. Instrument
- 31. Instrument
- 32. Instrument
- 33. Instrument
- 34. Instrument
- 35. Instrument
- 36. Instrument
- 37. Instrument
- 38. Instrument
- 39. Instrument
- 40. Instrument
- 41. Instrument
- 42. Instrument
- 43. Instrument
- 44. Instrument
- 45. Instrument
- 46. Instrument
- 47. Instrument
- 48. Instrument
- 49. Instrument
- 50. Instrument
- 51. Instrument
- 52. Instrument
- 53. Instrument
- 54. Instrument
- 55. Instrument
- 56. Instrument
- 57. Instrument
- 58. Instrument
- 59. Instrument
- 60. Instrument
- 61. Instrument
- 62. Instrument
- 63. Instrument
- 64. Instrument
- 65. Instrument
- 66. Instrument
- 67. Instrument
- 68. Instrument
- 69. Instrument
- 70. Instrument
- 71. Instrument
- 72. Instrument
- 73. Instrument
- 74. Instrument
- 75. Instrument
- 76. Instrument
- 77. Instrument
- 78. Instrument
- 79. Instrument
- 80. Instrument
- 81. Instrument
- 82. Instrument
- 83. Instrument
- 84. Instrument
- 85. Instrument
- 86. Instrument
- 87. Instrument
- 88. Instrument
- 89. Instrument
- 90. Instrument
- 91. Instrument
- 92. Instrument
- 93. Instrument
- 94. Instrument
- 95. Instrument
- 96. Instrument
- 97. Instrument
- 98. Instrument
- 99. Instrument
- 100. Instrument

Accessories and Application Aids

PROSEPT® Burs Container

Small soaking tray for the disinfection of instruments such as rotary instruments.



PROSEPT® Soaking Tray

Soaking tray for the disinfection and cleaning of medical instruments.



Trigger Sprayer I1

European-made trigger sprayer with 1.3 ml output for bottles with DIN 28/410 thread.



Canister Tap

Canister tap for 5 and 10 litre canisters with DIN 51 and DIN 61 thread.



Do you have any infection control related questions?
Let us know at www.udm.swiss

PROSEPT® Wipes Tub T1

Suitable for: PROSEPT® Wipes (14.5 x 20 cm, 120 pcs)



PROSEPT® Wipes Tub T2

Suitable for: PROSEPT® Wipes Jumbo (14.5 x 40 cm, 100 pcs),
PROSEPT® Wipes Value (14.5 x 20 cm, 250 pcs)



PROSEPT® Wipes Tub T3

Suitable for: PROSEPT® Wipes Premium (19.3 x 21 cm, 120 pcs)



Do you have any infection control related questions?
Let us know at www.udm.swiss



SWISS  MADE

UDM UNITED
DISINFECTANT
MANUFACTURERS

United Disinfectant Manufacturer AG
Allmendstrasse 21 • 8320 Fehraltorf • Switzerland
Tel +41 (0)55 511 03 33 • info@udm.swiss • www.udm.swiss