

# PROSEPT® Floor – Concentrated liquid for the disinfection and cleaning of non-invasive medical devices.



## PROSEPT® Floor

PROSEPT® Floor is a gentle but powerful concentrate for the disinfection and cleaning of large surfaces of non-invasive medical devices made of sensitive materials such as radiology equipments. The unique combination of active agents based on quaternary ammonium compounds result in a broad effective spectrum. PROSEPT® Floor is compatible with a wide range of materials such as acrylic glass, polycarbonate (PC) and PVC, does not leave sticky residues and is free of aldehydes, and phenols.

### Your Benefits

- ▶ **Extremely economical concentrate**  
The extremely economical concentrate significantly reduces your disinfection expenses.
- ▶ **Good material compatibility**  
The product is suitable for the disinfection of a wide variety of water-washable surfaces.
- ▶ **Simple application**  
Volumes are indicated on the lid for fast and easy dosing of disinfectant concentrates.

### Effective Spectrum / Exposure Times

Bactericidal	2% in 1 minute
Antibiotic-resistant bacteria	2% in 1 minute
Yeasticidal	2% in 1 minute
Antibiotic-resistant yeasts	2% in 1 minute
Enveloped viruses	2% in 5 minutes
Bovine viral diarrhea virus (BVDV)	2% in 5 minutes
Coronavirus	2% in 5 minutes
Hepatitis B virus	2% in 5 minutes
Hepatitis C virus	2% in 5 minutes
Herpes simplex virus	2% in 5 minutes
Human immunodeficiency virus (HIV)	2% in 5 minutes
Influenza A virus	2% in 5 minutes
Vaccinia virus	2% in 5 minutes
Limited virucidal	2% in 5 minutes
Adenovirus	2% in 1 minute
Norovirus	2% in 5 minutes

Bactericidal (EN 13727:2012+A1:2013\*) 2% in 1 min; yeasticidal (EN 13624:2013\*) 2% in 1 min; enveloped viruses (EN 14476:2013+A1:2015\*) 2% in 5 min; limited virucidal (EN 14476:2013+A1:2015\*) 2% in 5 min

^ Tested under low organic load / \* Tested under high organic load

### Dosing

Concentration	Concentrate + Water	Time
2% (= 1.0x Cap =20 ml)	20 ml+980 ml	5 min.

### Delivery Forms

OD-011060	2 litre bottle with measuring cap
-----------	-----------------------------------

### Easy to Use

- ▶ For your protection, wear a laboratory gown, household gloves and safety goggles.
- ▶ Prepare the working solution according to the dosing table.
- ▶ Completely wipe or immerse the surface of the medical device to be disinfected and leave to act.
- ▶ Adhere to contact time.
- ▶ If necessary, dry the surface with a clean single-use cloth.
- ▶ Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area.
- ▶ A list of compatible and incompatible materials is available on our website.
- ▶ If you are uncertain as to the material composition of the medical device, contact the manufacturer before use.

### Limitations

- ▶ FOR PROFESSIONAL USE ONLY. Not for use on invasive medical devices. Do not mix with other products.
- ▶ Make sure that all visible soiling has been removed before disinfection.

### Compatibility

Acrylic glass  
Aluminium  
Polyethylene  
Polypropylene  
Polyvinyl chloride  
Silicone rubber

Compatible  
Compatible with restriction  
Compatible  
Compatible  
Compatible  
Compatible

### Relevant Physical and Chemical Properties

Colour:	Green
Physical state:	Clear, slightly viscous liquid
Scent:	Lemon
pH-value:	6.0 - 8.0
pH-value in aqueous solution:	6.0 - 8.0
Flammability:	The product is not flammable
Relative density (g/cm <sup>3</sup> ):	0.9

### Composition in 100 g

100 g contains 10 g quaternary ammonium compounds, auxiliaries and water.

### Shelf Life

Unopened container: 36 months

The shelf life of the product is subject to correct storage and transportation conditions and the use of original containers. Keep containers tightly closed.

### Warning

H315 Causes skin irritation. H410 Very toxic to aquatic life with long lasting effects. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P332 + P313 If skin irritation occurs: Get medical advice/attention. P362 Take off contaminated clothing. P391 Collect spillage. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

### Your Local Distributor

Hygiene360 AG  
Büelstrasse 17  
8330 Pfäffikon ZH  
Switzerland  
Telephone: +41 (0)55 511 22 33  
info@hygiene360.com  
www.hygiene360.com

