

PROSEPT® Sensi – Alcohol-free, ready-to-use solution for the quick and comprehensive disinfection of non-invasive medical devices.



PROSEPT® Sensi

PROSEPT® Sensi is an alcohol-free, ready-to-use disinfection and cleaning product for non-invasive medical devices such as examination tables, dental treatment units, and medical appliances. The gentle formulation of PROSEPT® Sensi is particularly well suited for alcohol-sensitive surfaces made of PVC, silicone rubber, acrylic glass, or artificial leather. PROSEPT® Sensi offers a broad effective spectrum that includes pathogens such as MRSA and blood-borne viruses, including HBV, HCV, and HIV. PROSEPT® Sensi is free of volatile organic compounds (VOC), perfumes, and aldehydes, which are considered harmful to health.

Your Benefits

► Alcohol-free

The alcohol-free disinfectant is suitable for the disinfection of sensitive surfaces and objects such as acrylic glass.

► Long-lasting antimicrobial action

The antimicrobial action remains after the product has dried, thus, effectively reducing the risk of recontamination.

Effective Spectrum / Exposure Times

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 minutes
Mycobactericidal	60 minutes
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhoea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 minutes
Adenovirus	30 minutes
Norovirus	60 minutes
Poliovirus	60 minutes
Rotavirus	30 seconds

Bactericidal (EN 13727:2012+A2:2015*, EN 16615:2015*) 60 sec; tuberculocidal (EN 14348:2005*) 60 min; mycobactericidal (EN 14348:2005*) 60 min; yeasticidal (EN 13624:2013*, EN 16615:2015*) 60 sec; enveloped viruses (EN 14476:2013+A1:2015*) 15 sec; virucidal (EN 14476:2013+A1:2015*) 60 min

^ Tested under low organic load / * Tested under high organic load

Delivery Forms

OD-011004	200 ml foamer bottle
OD-011100	1 litre bottle
OD-011104	5 litre canister

Easy to Use

- ▶ For your protection, wear disposable gloves and safety goggles.
- ▶ Completely moisten and wipe the entire surface of the medical device.
- ▶ Leave disinfectant to act and adhere to contact time.
- ▶ Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area.
- ▶ A list of compatible and incompatible materials is available on our website.
- ▶ If you are uncertain as to the material composition of the medical device, contact the manufacturer before use.
- ▶ Wipe medical devices, which come into contact with skin, once a day with a water-soaked single-use cloth to prevent the build-up of biocides.

Limitations

- ▶ FOR PROFESSIONAL USE ONLY.
- ▶ Not for use on invasive medical devices.
- ▶ Do not mix with other products.
- ▶ Make sure that all visible soiling has been removed before disinfection.

Compatibility

Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Artificial leather	Compatible
Ceramic	Compatible
Chlorinated polyvinyl chloride	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Rhodium	Compatible
Silicone	Compatible
Stainless steel	Compatible
Steel	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

Relevant Physical and Chemical Properties

Colour:	Colourless
Physical state:	Clear, non-viscous liquid
Scent:	Apple
pH-value:	8.2 - 8.8
pH-value in aqueous solution:	Not applicable (ready-to-use solution)
Flammability:	The product is not flammable
Relative density (g/cm ³):	1.00

Composition

100 g contains 0.4 g benzalkonium chloride, 0.4 g N-alkyl aminopropyl glycine, auxiliaries and water.

Shelf Life

Unopened container:	48 months
---------------------	-----------

The shelf life of the product is subject to correct storage and transportation conditions and the use of original containers. Keep containers tightly closed.

No signal word according to Regulation (EC) 1272/2008

H412 Harmful to aquatic life with long lasting effects. P273 Avoid release to the environment. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/ container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

Your Local Distributor

United Disinfectant Manufacturers AG
 Allmendstrasse 21
 8320 Fehraltorf
 Switzerland
 Telephone: +41 (0)55 511 03 33
 info@udm.swiss
 www.udm.swiss

