



PROSEPT® Sonic – Liquid concentrate for the thorough cleaning of contaminated dental and surgical instruments in the ultrasonic cleaning device.



Delivery Forms

OD-061050

1 litre bottle

PROSEPT® Sonic

PROSEPT® Sonic is a non-foaming concentrate, based on high-performance alkali active components, for the efficient cleaning of dental and surgical instruments such as abrasors, amalgam carriers, spatulas, dental elevators, extraction forceps, and other precision instruments. PROSEPT® Sonic is free of aldehydes and phenols. PROSEPT® Sonic can be used for immersion procedures without the ultrasonic cleaning bath.

Your Benefits

- ▶ **Powerful cleanser**
The product's powerful cleansing properties make pre- and post-cleaning unnecessary.
- ▶ **Excellent cleansing properties**
The disinfectant leaves no streaks and therefore makes additional wiping redundant.
- ▶ **Simple application**
The bottle cap contains volume markings for fast and easy dosing.



Easy to Use

- ▶ For your protection, wear a laboratory gown, household gloves and safety goggles.
- ▶ Prepare the working solution according to the dosing table.
- ▶ Use an ultrasonic cleaning device for 30 minutes.
- ▶ After cleaning, rinse instruments thoroughly under running water of at least drinking quality for 10 seconds.
- ▶ Replace working solution daily or when visibly contaminated.
- ▶ Adhere to reprocessing recommendations of the instrument manufacturer.
- ▶ Not suitable for the cleaning of color-coded instruments and instruments that are fully or partly made of natural bristles, nylon, or aluminium.
- ▶ A list of compatible and incompatible materials is available on our website.
- ▶ If you are uncertain as to the material composition of the medical device, contact the manufacturer before use.

Limitations

- ▶ FOR PROFESSIONAL USE ONLY.
- ▶ Do not mix with other products.

Compatibility

Acrylic glass
Acrylonitrile butadiene styrene
Chlorinated polyvinyl chloride
Polybutylene terephthalate
Polycarbonate
Polyethylene
Polyethylene terephthalate
Polyphenylene oxides
Polyphenylsulfone
Polypropylene
Polystyrene
Polysulfone
Polyvinyl chloride
Polyvinylidene fluoride
Rhodium
Silicone
Thermoplastic elastomer
Thermoplastic polyurethane
Ultra-high-molecular-weight-polyethylene

Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible
Compatible

Relevant Physical and Chemical Properties

Colour:	Light yellow
Physical state:	Clear, non-viscous liquid
Scent:	Aromatic
pH-value:	14.0
pH-value in aqueous solution:	12.5 - 13.5
Flammability:	The product is not flammable
Relative density (g/cm ³):	1.2

Composition

100 g contains 16 g potassium hydroxide, corrosion inhibitors, auxiliaries and water.

Shelf Life

Unopened container:	60 months
---------------------	-----------

The shelf life of the product is subject to correct storage and transportation conditions and the use of original containers. Keep containers tightly closed. The shelf life of a working solution in use depends on the degree of contamination. In case of heavy use, a daily change is recommended. Always replace the working solution if it becomes turbid or changes colour.

Danger

H290 May be corrosive to metals. H314 Causes severe skin burns and eye damage. P260 Do not breathe fume/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

Your Local Distributor

United Disinfectant Manufacturers AG
Allmendstrasse 21
8320 Fehraltorf
Switzerland
Telephone: +41 (0)55 511 03 33
info@udm.swiss
www.udm.swiss

