



## PROSEPT® Spray – Ready-to-use solution for the quick and comprehensive disinfection of non-invasive medical devices.



### PROSEPT® Spray

PROSEPT® Spray is a ready-to-use disinfection and cleaning product for non-invasive medical devices such as examination tables, dental treatment units, and medical appliances. Thanks to its time-proven combination of ethanol and 1-propanol, PROSEPT® Spray offers a broad effective spectrum and a very short contact time. The low content of volatile organic compounds (VOC) of less than 50% prevents cracking of the disinfected surfaces and reduces skin and breathing problems common with disinfectants containing a high alcohol percentage. PROSEPT® Spray is free of aldehydes and phenols and is fragranced with a pleasant lemon perfume.

### Your Benefits

- ▶ **Alcohol content of less than 50%**  
The disinfectant is particularly suitable for use on sensitive surfaces such as artificial leather covers.
- ▶ **Ergonomically designed bottle**  
The ergonomically designed bottle ensures a safe and easy application of the disinfectant.

### Effective Spectrum / Exposure Times

Bactericidal	60 seconds
Antibiotic-resistant bacteria	60 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	60 seconds
Antibiotic-resistant yeasts	60 seconds
Virucidal against enveloped viruses	15 seconds
Bovine viral diarrhoea virus (BVDV)	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
Human immunodeficiency virus (HIV)	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Virucidal	60 seconds
Adenovirus	60 seconds
Norovirus	30 seconds
Poliovirus	60 seconds

Bactericidal (EN 13727:2012+A2:2015\*, EN 16615:2015\*) 60 sec; tuberculocidal (EN 14348:2005\*) 60 sec; mycobactericidal (EN 14348:2005\*) 60 sec; yeasticidal (EN 13624:2013\*, EN 16615:2015\*) 60 sec; enveloped viruses (EN 14476:2013+A1:2015\*) 15 sec; virucidal (EN 14476:2013+A1:2015\*) 60 sec

^ Tested under low organic load / \* Tested under high organic load

### Delivery Forms

OD-041005	250 ml bottle
OD-041013	1 litre bottle
OD-041025	5 litre canister

### Easy to Use

- ▶ For your protection, wear disposable gloves and safety goggles.
- ▶ Completely moisten and wipe the entire surface of the medical device.
- ▶ Leave disinfectant to act and adhere to contact time.
- ▶ Use on alcohol compatible materials only.
- ▶ Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area.
- ▶ A list of compatible and incompatible materials is available on our website.
- ▶ If you are uncertain as to the material composition of the medical device, contact the manufacturer before use.

### Limitations

- ▶ FOR PROFESSIONAL USE ONLY.
- ▶ Not for use on invasive medical devices.
- ▶ Do not mix with other products.
- ▶ Make sure that all visible soiling has been removed before disinfection.

### Compatibility

Acetal polymer	Compatible
Acrylic glass	Not compatible
Acrylonitrile butadiene styrene	Compatible
Aluminium	Compatible
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Rhodium	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

### Relevant Physical and Chemical Properties

Colour:	Colourless
Physical state:	Clear, non-viscous liquid
Scent:	Of alcohol, lemon
pH-value:	Neutral
pH-value in aqueous solution:	Not applicable (ready-to-use solution)
Flammability:	The product is flammable
Relative density (g/cm <sup>3</sup> ):	0.92

### Composition

100 g contains 20 g ethanol, 28 g 1-propanol, 0.056 g quaternary ammonium compounds, auxiliaries and water.

### Shelf Life

Unopened container:	48 months
---------------------	-----------

The shelf life of the product is subject to correct storage and transportation conditions and the use of original containers. Keep containers tightly closed.

### Danger

H226 Flammable liquid and vapour. H318 Causes serious eye damage. H336 May cause drowsiness or dizziness. P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P261 Avoid breathing fume/mist/vapours/spray. P280 Wear protective gloves/protective clothing/eye protection/face protection. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

### Your Local Distributor

United Disinfectant Manufacturers AG  
 Allmendstrasse 21  
 8320 Fehraltorf  
 Switzerland  
 Telephone: +41 (0)55 511 03 33  
 info@udm.swiss  
 www.udm.swiss

