

## PROSEPT® Floor – Concentrated liquid for the disinfection and cleaning of non-invasive medical devices.



### PROSEPT® Floor

PROSEPT® Floor is a gentle but powerful concentrate for the disinfection and cleaning of large surfaces of non-invasive medical devices made of sensitive materials such as radiology equipment. The unique combination of active agents based on quaternary ammonium compounds result in a broad effective spectrum. PROSEPT® Floor is compatible with a wide range of materials such as acrylic glass, polycarbonate (PC) and PVC, does not leave sticky residues and is free of aldehydes and phenols.

### Your Benefits

- ▶ **Extremely economical concentrate**  
The extremely economical concentrate significantly reduces your disinfection expenses.
- ▶ **Good material compatibility**  
The product is suitable for the disinfection of a wide variety of water-washable surfaces.
- ▶ **Simple application**  
Volumes are indicated on the lid for fast and easy dosing of disinfectant concentrates.

### Effective Spectrum / Exposure Times

Bactericidal	2% in 1 minute
Antibiotic-resistant bacteria	2% in 1 minute
Yeasticidal	2% in 1 minute
Antibiotic-resistant yeasts	2% in 1 minute
Virucidal against enveloped viruses	2% in 5 minutes
Bovine viral diarrhoea virus (BVDV)	2% in 5 minutes
Coronavirus	2% in 5 minutes
Hepatitis B virus	2% in 5 minutes
Hepatitis C virus	2% in 5 minutes
Herpes simplex virus	2% in 5 minutes
Human immunodeficiency virus (HIV)	2% in 5 minutes
Influenza A virus	2% in 5 minutes
Vaccinia virus	2% in 5 minutes
Limited spectrum virucidal	2% in 5 minutes
Adenovirus	2% in 1 minute
Norovirus	2% in 5 minutes

Bactericidal (EN 13727:2012+A1:2013\*, EN 16615:2015^\*) 2% in 1 min; yeasticidal (EN 13624:2013\*, EN 16615:2015^\*) 2% in 1 min; enveloped viruses (EN 14476:2013+A1:2015\*) 2% in 5 min; limited virucidal (EN 14476:2013+A1:2015^\*) 2% in 5 min  
^ Tested under low organic load / \* Tested under high organic load

### Dosing

Concentration	Concentrate + Water	Time
2% (= 1.0x Cap = 20 ml)	20 ml + 980 ml	5 min.

### Delivery Forms

OD-011060	2 litre bottle
OD-010065	5 litre canister

### Easy to Use

- ▶ For your protection, wear a laboratory gown, household gloves and safety goggles.
- ▶ Prepare the working solution according to the dosing table.
- ▶ Completely wipe or immerse the surface of the medical device to be disinfected and leave to act.
- ▶ Adhere to contact time.
- ▶ If necessary, dry the surface with a clean single-use cloth.
- ▶ Before first time use, test material compatibility by applying the disinfectant to a small inconspicuous area.
- ▶ A list of compatible and incompatible materials is available on our website.
- ▶ If you are uncertain as to the material composition of the medical device, contact the manufacturer before use.

### Limitations

- ▶ FOR PROFESSIONAL USE ONLY. Not for use on invasive medical devices.
- ▶ Do not mix with other products.
- ▶ Make sure that all visible soiling has been removed before disinfection.

### Compatibility

Acetal polymer	Compatible
Acrylic glass	Compatible
Acrylonitrile butadiene styrene	Compatible
Aluminium	Compatible with restriction
Artificial leather	Compatible
Chlorinated polyvinyl chloride	Compatible
Polyamide	Compatible
Polybutylene terephthalate	Compatible
Polycarbonate	Compatible
Polyetherimide	Compatible
Polyethylene	Compatible
Polyethylene terephthalate	Compatible
Polyphenylene oxides	Compatible
Polyphenylsulfone	Compatible
Polypropylene	Compatible
Polystyrene	Compatible
Polysulfone	Compatible
Polyvinyl chloride	Compatible
Polyvinylidene fluoride	Compatible
Silicone	Compatible
Thermoplastic elastomer	Compatible
Thermoplastic polyurethane	Compatible
Ultra-high-molecular-weight polyethylene	Compatible

### Relevant Physical and Chemical Properties

Colour:	Green
Physical state:	Clear, slightly viscous liquid
Scent:	Lemon
pH-value:	6.0 - 8.0
pH-value in aqueous solution:	6.0 - 8.0
Flammability:	The product is not flammable
Relative density (g/cm <sup>3</sup> ):	0.99

### Composition

100 g contains 10 g quaternary ammonium compounds, auxiliaries and water.

### Shelf Life

Unopened container:	36 months
---------------------	-----------

The shelf life of the product is subject to correct storage and transportation conditions and the use of original containers. Keep containers tightly closed.

### Warning

H315 Causes skin irritation. H410 Very toxic to aquatic life with long lasting effects. P280 Wear protective gloves/protective clothing/eye protection/face protection. P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P332 + P313 If skin irritation occurs: Get medical advice/attention. P362 Take off contaminated clothing. P391 Collect spillage. P501 Dispose of contents/container in accordance with local and national regulations. For more information, please refer to Safety Data Sheet available on our website.

### Your Local Distributor

United Disinfectant Manufacturers AG  
 Allmendstrasse 21  
 8320 Fehraltorf  
 Switzerland  
 Telephone: +41 (0)55 511 03 33  
 info@udm.swiss  
 www.udm.swiss

